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TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY. — Nebezpečenstvo udusenía. Nepoužívajte toto

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Symbol  znamená UPOZORNENIE.

INSTRUCTIONS FOR USE

Dear customer, thank you for purchasing our product. Please read the operating instructions carefully before putting the appliance into operation and keep these instructions including the warranty, the receipt and, if possible, the box with the internal packing.

I. GENERAL SAFETY INSTRUCTIONS

- Consider the instructions for use as a part of the appliance and pass them on to any other user of the appliance.
- Children and persons with limited physical or mental capabilities or lack of experience may only use this appliance if they are supervised or have been instructed of the use of the appliance in a safe manner and understand possible risks. Unattended children must not perform cleaning and maintenance intended to be performed by users. Children must not play with the appliance.
- Never use the appliance if it does not correctly work or fell on the ground and is damaged.
- Use zinc-carbon or alkaline battery of AAA size. Observe the indicated polarity when exchanging batteries. If you are going not to use the toothbrush for a long time, remove the batteries.
- If the battery shows any indications of corrosion or leaking electrolyte, immediately replace it, otherwise the brush might get damaged.
- Dispose the empty battery in the appropriate manner (see article **V. ECOLOGY**).
- Keep the battery and toothbrush out of reach of children and incapacitated persons. Person who has swallowed the battery must immediately seek for medical assistance!
- The appliance is only intended for individual use at home for personal oral cavity care! It is not intended for commercial use, for more than one patient in dentist's consulting rooms or in any other healthcare facilities!
- Do not use the application in exterior or where aerosol products are used (such as **sprays**), or where oxygen is handled.
- **Do not immerse the driver unit in water or other liquids!**
- Do not insert the appliance or its accessories in any body cavities, except for oral cavity.
- The appliance is not designed to care for animal teeth.
- Clean the accessories after each use. Failure to adhere to this instruction might cause worsening your state of health, because without disinfection, spread of mildew and bacteria is supported.
- **Do not leave the appliance unattended!**
- Protect the appliance from dust, chemicals, toiletries, liquid cosmetics, large fluctuations in temperature, and direct sunshine.
- Do not attempt to repair the appliance by yourselves.
- Keep the accessories out of reach of children and incapacitated persons.

SPECIAL SAFETY INSTRUCTIONS

- If your toothpaste contains peroxide, baking soda or other bicarbonate (frequently used in whitening toothpastes), thoroughly clean the brush and the handle after use to prevent damaging the appliance.
- Do not use the appliance in care for children, incapacitated persons or persons who are unable to express their feelings.
- If you or your child have undergone tooth or gum operation recently, suffered long-term problems with teeth or were subject to extensive medical interventions in the last 2 months, or if you suspect problems of dental character, consult paediatric dentist before using the appliance.
- If you experience any health troubles (for example considerable bleeding, or if any bleeding persists even after 2 weeks of use of the brush), stop using the appliance and contact your physician for consultation.
- The appliance complies with the safety standards for electromagnetic radiation. However if you or your child uses any of the following devices: cardiostimulator, electrocardiograph or other medical device implanted in your body or used as a support, consult your physician or the implant manufacturer before use.
- For hygienic reasons do not share the toothbrush with another person or more persons. Infections could spread or oral cavity inflammation could be initiated.
- Only use the appliance with the original accessories designated for the respective type. Do not use damaged/deformed adapters.
- Only use this appliance including its accessories for the designated purpose as described in this instruction manual. The appliance is only designed for cleaning teeth, gums and tongue. Never use the appliance for any other purpose.
- **WARNING:** There is a risk of injury in case of incorrect use (not in accordance with the manual) of the appliance.
- All text in other languages and pictures shown on the packaging, or the product itself, are translated and explained at the end of this language mutation.
- The producer does not undertake any responsibility for damage caused by improper use of the appliance (e.g. **injury etc.**) and its guarantee does not cover the appliance in the case of failure to comply with the safety warning above.

II. DESCRIPTION OF FEATURES

A – Brush body

A1 – ON/OFF switch

A2 – battery compartment

B – Brush (2 pcs)

B1 – brush head

B2 – coloured ring to distinguish users

C – Battery cover

Note: Replacement brushes can be purchased under order number ETA071190100.

III. INSTRUCTIONS FOR USE

Remove all packaging materials and take out the appliance and any accessories. Parts which come in contact with the oral cavity must be rinsed in warm drinking water before first use.

INSERTING THE BATTERY

Remove the battery cover **C** as per Fig. 1 and insert a new AAA battery (included) into the compartment **A2**. Close the cover in the opposite way. Depending on the cleaning mode and on the battery used, the brush will provide approximately 3-6 hours of operation which corresponds to about 63 months when used twice a day for two minutes.

USE

Before use, connect the brush **B** with the brush body **A** as per Fig. 1. To remove the brush, proceed the other way round. **Do not switch on the appliance without the brush attached!** Moisten the brush and apply a small amount of toothpaste on it. Put the toothbrush to the child's teeth at a slight angle towards the gums (see Fig. 2A) and switch it on by pressing the switch button **A1**. When the appliance is switched on, move the brush slowly along the teeth in small circular motions, applying a slight pressure, so as not to miss any tooth or area in the mouth. Clean each tooth for about 2-3 seconds. Each part of the mouth (quadrant) should be cleaned for about 30 seconds (see Fig. 2B), a short interruption of vibrations on the brush will prompt you to move to another quadrant after this time elapses. Switch off the appliance after use by pressing button **A1**.

The brush operating time is limited to 2 minutes which corresponds to the time for cleaning recommended by experts. If you wish to continue cleaning, simply switch on the brush again.

Notes

- If this is the first sonic toothbrush your child is using and the child has not had any experience with this cleaning technology, it might take them a while to get used to the specific vibrations (tingling). However, most users say that they got used to it quickly after a few uses and cleaning became perfectly natural for them.
- Do not apply excessive pressure on your teeth.
- You can also use this brush when you have fixed braces, facets, crowns, etc. However, please bear in mind that if the brush is used with braces, the brush bristles are subject to more wear.
- We recommend that you replace the brush **B** with a new one after 3 months of use. However, if the brush is showing signs of significant wear or damage, replace it sooner.

IV. CLEANING AND MAINTENANCE

For maximum hygiene, clean the brush B after each use! Do not use isopropyl rubbing alcohol, vinegar or bleach to clean brush body as this may cause discoloration. Rinse the brush head **B1** and the brush body **A** under running water and wipe the brush body with a damp cloth. Do not immerse the brush body **A** in water. Do not use a dishwasher for cleaning! Allow the brush body **A** to dry properly before its re-use.

Store the appliance and all the accessories in a dry, dust-free and safe place, out of reach of children and incapacitated persons.

V. ENVIRONMENTAL PROTECTION

If dimensions permit, all pieces have printed symbols of materials used for production of packing, components and accessories as well as recycling symbols.

The symbols on the product or documentation means the electric or electronic parts used in the appliance may not be disposed off together with communal waste. Hand over the electric or electronic part to a waste collection yard for free disposal. Proper disposal of the product will help preserving valuable natural resources and preventing potential negative impacts on environment and human health that would otherwise result from improper waste disposal. For more details, please contact your local authorities or collection yard. Fines may be imposed for improper disposal of this type of waste pursuant to national directives. Remove battery from the appliance and dispose it off in a suitable way using special collection centres. Do not incinerate the battery!

VI. TECHNICAL DATA

Weight (g) approximately	30
Protection class of the appliance	III.
Size, (mm):	170 x 25 x 25

The manufacturer reserves the right to change the technical specifications and accessories for the respective models. Expert service must do the large maintenance or to intervene into inside parts of the appliance! Non-observance of the instructions of the manufacturer cancels the right for guarantee repair!

WARNINGS AND SYMBOLS USED ON THE APPLIANCE, PACKAGING OR IN THE INSTRUCTIONS MANUAL:

HOUSEHOLD USE ONLY; DO NOT IMMERSE IN WATER OR OTHER LIQUIDS;

TO AVOID DANGER OF SUFFOCATION, KEEP THIS PLASTIC BAG AWAY FROM BABIES AND CHILDREN. DO NOT USE THIS BAG IN CRIBS, BEDS, CARRIAGES OR PLAYPENS. THIS BAG IS NOT A TOY.

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© DATE 18/4/2018